



JBone-03K

AI based Bone Age Analysis Solution

| Summary

Bone age analysis is required for prediction of the precocious puberty, the growth rate & adult height.

JBone-03K solution assists in diagnosing precocious puberty or growth disorders by measuring bone age within 8 months of error without any other clinical information.

| Key Components & Performance

- X-ray image upload
- X-ray image visualization
- Area for Bone age calculation auto-detection and visualization
- Bone age auto-calculation

| Input / Output

- **Input Data** Patient's X-Ray DICOM
- **Output Information** Bone Age analysis information